



LEADERS IN SAMPLE INTRODUCTION TECHNOLOGY

5 0 0 0



solids pyrolyzer

pulse pyrolysis

programmed pyrolysis

reactant gas

multi-step pyrolysis

Pyroprobe® 5000

The CDS Pyroprobe® Model 5000



is the most advanced thermal sample preparation instrument available. Platinum filaments are rapidly heated for pulse pyrolysis work, or slowly heated with controlled rates for programmed analyses. Set temperatures in 1°C increments to 1400°C for a wide heating range and more precise pyrolysis temperatures. Analytical runs may be programmed for up to five steps per sample, with automatic control of the on-line valve, interface temperature, GC ready sensing and GC start for each step. When configured

with the trapping option, the Pyroprobe® may be used to collect analytes from slow rate pyrolysis, thermal desorption or reactant gas pyrolysis. The unique interfacing design permits a direct pyrolysis path to the GC inlet or rapid sample heating and transfer to the trap without interrupting the pneumatics of the GC. Simple PC control of all parameters is presented in an easy to program window, permitting method development, storage, notation and editing.

Let our team of scientists help you with your analytical challenges

automatic

product features



specify up to five runs per sample—
automatically

Window-based CDS 5000 DCI Control

platinum filament pyrolyzer for pulse-pyrolysis and multistep, programmed pyrolysis

simplified sample loading

programmable interface

easy transfer line connection to GC

PC control

three modes of operation

run, dry, clean—all user selectable

built-in trapping zone for:

thermal desorption

reactant gas operation

slow rate pyrolysis with trapping

interface directly to GC or to trap

product specifications

Model 5000

Pulse Pyrolysis

Filament Temperature: Programmable in 1°C increments to 1400°C
Heating Rates: 0.01 to 20.0°C/ms (10 to 20,000°C/second)
Clean and Dry: User selectable

Programmed Pyrolysis

Heating Rates: 0.01 to 999.9°C/sec and 0.01 to 999.9°C/min
Multiple Steps: Each step may include GC start

Interface

1500 Valved Interface: Isothermal temperature, set in 1°C increments



control

advanced capabilities

Model 5150

all above, plus

Multiple steps: Each interface step includes initial, ramp and final
Each step may include Pyroprobe® and GC start

Interface: Low-mass programmable zone on 5150 unit
Temperature: Settable in 1°C/min increments to 350°C
Heating Rates: Programmable in 1°C/min increments to 60°C/min


Heated Zones: 4 auxiliary zones
Valve Oven: Settable in 1°C increments to 300°C
Transfer Line: Settable in 1°C increments to 325°C



Model 5200

for headspace trapping and reactant gas pyrolysis

Interface: Pyrolysis to built-in trap or direct to GC
Temperature: Settable in 1°C increments to 350°C
Heating Rates: Programmable in 1°C/min increments to 60°C/min
Trapping Tube: Ambient to 350°C



Overall Dimensions: 23 cm W x 24 cm H x 31 cm D
(5000 without interface)





LEADERS IN SAMPLE INTRODUCTION TECHNOLOGY


MAKERS OF CDS MODEL 8000 SAMPLE CONCENTRATOR AND DYNATHERM THERMAL DESORPTION SYSTEMS

CDS Analytical, Inc. The World's Most Versatile Pyrolysis, Headspace, and Thermal Desorption Systems

About CDS	Products	Application Notes	Customer Support	Sales Distribution	Links	What's Happening	Quote/ Lit. Contact	Contact Us
---------------------------	--------------------------	-----------------------------------	----------------------------------	------------------------------------	-----------------------	----------------------------------	-------------------------------------	----------------------------

Welcome to CDS Analytical
 For over 30 years, CDS Analytical, Inc. has been providing analytical instruments worldwide for the analysis of complex samples by GC, GC/MS, and FT-IR. Our products include pyrolyzers, thermal desorbers, autosamplers and sample concentrators for pyrolysis, thermal desorption, dynamic headspace, and trace analysis. We pride ourselves on providing application specific solutions, backed by the best service and support in the industry.

Dynatherm Product Line  **Application Lab Support** 



465 Limestone Rd
 Oxford, PA 19363
 Tel: +1.610.932.3636
 Fax: +1.610.932.4158

Pharmaceutical Polymers Forensics Electronics
 Environmental Paper Homeland Other

for more information

on CDS products and

applications, go to

www.cdsanalytical.com

Worldwide Headquarters
465 Limestone Road
P.O. Box 277
Oxford, PA 19363-0277 USA
Toll Free 800.541.6593
610.932.3636
Fax 610.932.4158