


## Thermal Desorption of Gasoline Range Hydrocarbons from Soil using a Pyroprobe


Soil may be screened easily for contamination by placing a small sample into a quartz tube and using the Pyroprobe coil filament to desorb the organics from the soil matrix. Care must be taken not to heat the soil to a temperature high enough that the natural organic content is pyrolyzed and produces volatiles that complicate the analysis. In this example, it was determined that at 300°C, the soil organics produced volatiles, including furans and phenolics, but at 200°, the blank was relatively clean.

Figure 1 shows the result of heating 20 mg of dry soil to 200°C, after hydrocarbons had been added to make the total concentration 100 PPM. Sample preparation consisted only of weighing the sample and then transferring it to a quartz tube using a mini-vacuum pump. The soil was placed in the center half of the tube, and held in place using plugs of quartz wool. The quartz tube was then placed into the coil of the Pyroprobe, and heated in the interface of the Pyroprobe 5200. The resulting volatilized hydrocarbons were collected onto the sorbent trap before desorption to the GC.

Figure 2 shows a blank run of the same soil without the added hydrocarbons, at 200°C. A few small peaks are seen at this temperature, which also are present in the run of the contaminated soil, but they do not mask or interfere with the identification of the hydrocarbons resulting from the contamination.



**Figure 1.** Hydrocarbons in soil - 100 PPM total.



**Figure 2.** Blank soil at 200°C.

## Equipment

### Pyroprobe 5200

Interface:	200°C for 4 minutes
Filament:	200°C for 30 seconds
Valve Oven:	325°C
Transfer line:	325°C
Trap desorption:	325°C for 4 minutes
Dry:	350°C for 60 seconds
Clean:	800°C for 10 seconds

### GC/MS

Column:	30 m x 0.25 mm 5% phenyl methyl silicone
Carrier:	Helium
Split:	50:1
Oven program:	40°C for 2 minutes 10°C/min to 300 °C

FOR MORE INFORMATION  
CONCERNING THIS APPLICATION,  
WE RECOMMEND THE  
FOLLOWING READING:

D. White and L. Beyer, J. Anal. Appl.  
Pyrolysis, 50 (1999) 63-76.

© courtesy of CDS Analytical Inc.

### Headquarters

JSB International  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 251 47 53  
F +31 (0) 40 251 47 58

Zoex Europe  
Tramstraat 15  
5611 CM Eindhoven  
T +31 (0) 40 257 39 72  
F +31 (0) 40 251 47 58

### Sales and Service

Netherlands  
Apolloweg 2B  
8239 DA Lelystad  
T +31 (0) 320 87 00 18  
F +31 (0) 320 87 00 19

Belgium  
Grensstraat 7  
Box 3 1831 Diegem  
T +32 (0) 2 721 92 11  
F +32 (0) 2 720 76 22

Germany  
Max-Planck-Strasse 4  
D-47475 Kamp-Lintfort  
T +49 (0) 28 42 9280 799  
F +49 (0) 28 42 9732 638

UK & Ireland  
Cedar Court,  
Grove Park Business Est.  
White Waltham, Maidenhead  
Berks, SL6 3LW  
T +44 (0) 16 288 220 48  
F +44 (0) 70 394 006 78

[info@go-jsb.com](mailto:info@go-jsb.com)  
[www.go-jsb.com](http://www.go-jsb.com)

