

Agilent 7820A Gas Chromatograph System

Data Sheet

The Agilent 7820A gas chromatograph, inheriting Agilent's legendary expertise in GC and proven quality as industry leader, generates reliable results with the minimized complexity for customers' routine analyses, run after run, day after day.

The simplified GC front panel keys and display provide sequence information, instrument conditions, and run status. Full electronic pneumatics control (EPC) is available for all inlets and detectors.

Configurable with either a 50-vial injector (as option) or a 16-vial injector.

A flexible software choice is provided. From OpenLAB CDS ChemStation Edition or EZChrom Edition, OpenLAB CDS ChemStation VL or EZChrom VL, to OpenLAB CDS EZChrom Compact, everything is designed to help you make the most of every run, and every workday.

Dimensions and Weight

 Height
 49 cm

 Width
 56 cm

 Depth
 51 cm

 Average weight
 50 kg

GC Front Panel Keys and Display

Available in English or Chinese

Environmental Conditions

Indoor use

Ambient operating temperature 15 to 30 °CAmbient operating humidity 30 to 70%Storage extremes -40 to 70 °COperating altitude 3,100 m



Safety and Regulatory Certifications

Safety Standards

Canadian Standards Association (CSA) C22.2 No. 61010
CSA/Nationally Recognized Test Laboratory (NRTL) UL61010
International Electrotechnical Commission (IEC) 61010
EuroNorm (EN) EN61010

Electromagnetic compatibility (EMC) and radio frequency interference (RFI) regulation conformity

CISPR 11/EN 55011 Group 1, Class A

IEC/EN 61326

Designed and manufactured under a quality system registered to ISO 9001. The Declaration of Conformity is available.

System Overall Performance*

* Using 7820A with EPC (splitless), ALS, and Agilent Data System for analysis of tridecane (2 ng to the column). Results may vary with other samples and conditions.

Retention time repeatability < 0.06% Peak area repeatability < 2%

Power Requirements

100 V (+10%, -10%) 120 V (+10%, -10%) 200 V (+10%, -10%) 220 V (+10%, -10%) 230 V (+10%, -10%) 240 V (+10%, -10%)

Frequency 47.5~63 Hz

1,500 W (max) at 100 V, 2,250 W (max) at all other voltages

Column Oven

Dimensions	$28.0 \times 30.5 \times 16.5 \text{ cm}$		
Operating temperature	8 °C above ambient to 425 °C		
Temperature setpoint resolution	1 °C		
Maximum temperature ramp rate	75 °C/min (see Table 1)		
Maximum run time	999.99 min		
Temperature programming ramps	5		
Ambient rejection	< 0.01 °C per 1 °C		
Oven temperature ramp	≤ 2%		
Programming temperature	< 1%		

Typical heating-up profile and cooldown rate are shown in Figures 1 and 2.

Heated Zones

- Five independent heated zones, not including oven (two inlets, two detectors, and one auxiliary)
- · 350 °C Maximum operating temperatures for auxiliary zone
- · Support up to two heated valves




Figure 1. Typical oven heat up profile.




Figure 2. Typical oven cooldown profile.

Table 1. Typical 7820A GC Oven Ramp Rates

Temperature (°C)	220 V oven rates (°C/min)	
50 to 70	75	
70 to 115	45	
115 to 175	40	
175 to 300	30	
300 to 425	20	

For 100 V oven, the maximum temperature is 350 °C with a maximum ramp rate of 30 °C/min.

Electronic Pneumatics Control (EPC)

Available on all inlet and all detectors

Inlet Module

Pressure setpoint and control

precision 0.01 psi or 0.069 KPa

 $< \pm 5\%$ Flow sensor accuracy

Detector Module Accuracy

8% of setpoint

Inlets

Maximum inlets installed two

Inlets available Purged packed (PP)

Split/splitless capillary (S/SL)

Purged Packed

· Electronic flow control

· 400 °C maximum operating temperature

Maximum flow < 100 mL/min

Adapters included for 1/4-inch and 1/8-inch packed columns and for

0.530-mm capillary columns

Electronic pressure/flow control

Maximum operating temperature 400 °C

0 to 60 psi or 0 to 413.69 KPa Pressure range

250:1 Maximum split ratio

Total flow setting range

0 to 200 mL/min $\rm N_2$ 0 to 500 mL/min $\rm H_2$ or He

Detectors

All detectors use electronic pneumatic control (EPC) for detector gases. Up to two detectors may be installed.

Available Detectors

· Flame ionization detector (FID)

· Thermal conductivity detector (TCD)

Electron capture detector (ECD)*

Nitrogen phosphorous detector (NPD)

· Flame Photometric detector (FPD)

FID

Electronic pressure/flow control

Maximum operating temperature 425 °C

< 3 pg carbon/s as tridecane MDL

 $> 10^7$ range with N₂ carrier and Linear dynamic range

0.29-mm id jet

Maximum data acquisition rate 100 Hz **TCD**

Electronic pressure/flow control

Maximum operating temperature 400 °C

< 800 pg propane/mL using He carrier (MDL may be affected by laboratory

environment)

Linear dynamic range 10⁵ (± 10%)

ECD*

Electronic pressure/flow control

Equipped with hidden anode and high-velocity flows for contamination resis-

tance

Maximum operating temperature 400 °C

Makeup gas types argon/5% methane or nitrogen

< 15 mCi 63Ni Radioactive source MDI < 0.02 pg/mL lindane > 104 with lindane Dynamic range

50 Hz maximum data acquisition rate

*ECD not supported in Japan

NPD

Electronic pressure/flow control

Maximum operating temperature 400 °C

MDL < 0.4 pg N/s, < 0.2 pg P/s with

azobenzene/malathion mixture

25,000 to 1 gN/gC, 75,000 to 1 gP/gC Selectivity

with azobenzene/malathion mixture

 $> 10^4 \text{ N}, > 10^4 \text{ P with}$ Dynamic range

azobenzene/malathion mixture

Data acquisition rate up to 100 Hz

FPD

Single wavelength

MDL < 200 fg P/s, < 6 pg S/s with

methylparathion

 $> 5 \times 10^2$ S, 10^4 P with Dynamic range

methylparathion

Selectivity 10^6 g S/g C, 10^6 g P/g C

Data acquisition rate up to 200 Hz

Standard EPC for three gases Air 0 to 200 mL/min 0 to 250 mL/min H.

> Makeup gas 0 to 130 mL/min

Maximum operating temperature 250 °C

Optional ALS

Supports one 7650A autoinjector with capacity for

50 sample vials

Supports one 7693A autoinjector with capacity for

16 sample vials

Data Communications

- · One analog output channel (1 mV, 1 V, and 10 V output available) as standard
- · Remote start/stop
- LAN

OpenLAB CDS Software Platform Hardware and Software Requirements

OpenLAB CDS software platform is necessary for complete 7820 system control, data analysis, and reporting.

The table below lists the hardware and software requirements for OpenLAB CDS software platform installations.

Product	Operating system	CPU (minimum)	RAM (minimum available to EZChrom)	Disk space (minimum available to EZChrom)	Ports available
OpenLAB CDS Workstation	Windows XP Professional (SP2) Internet Explorer 6.0 or 7.0 Microsoft.NET version 3.0	1.0 GHz	2.0 GB	10 GB	One Ethernet LAN port per system

- English version of OpenLAB CDS is validated on Chinese, Japanese, English, and Western European language versions of Windows.
- Japanese version of OpenLAB CDS is validated on Japanese and English language versions of Windows
- · Chinese version of OpenLAB CDS is validated on Chinese and English language versions of Windows.

GC Operation Without OpenLAB CDS Software

If an integrator or a third-party software package is being utilized for data analysis, a software virtual keyboard will be shipped with the 7820 GC. This standalone utility can run on a PC and allow creation of GC methods, which are then downloaded to the GC via a LAN cable.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2012 Published in the USA April 13, 2012 5990-3279EN

