

LC/MS Application Note







LIFE SCIENCE


Infusion Automation - Expert infusion automation for your high sensitivity metabolomics and lipidomics experiments

TRAJAN Automation

LEAP >>>

High-throughput, high sensitivity and low carryover: the new infusion automation from Trajan LEAP Automation offers you everything you need to bring your lipidomics mass spec workflow to its peak performance.

Move from manual sample infusion to the highest standard of infusion automation. With rapid cleanup between injections and best-in-class software architecture to integrate into your workflow, our infusion automation will elevate your lipidomics experiments.



An exclusive inert configuration is available for your trickiest lipid classes - ensuring your samples never come into contact with any polymers.

> "We just ran a whole pile of samples. It ran 24/7 without a hitch. Very reliable."

Calvin Vary, PhD Proteomics and Lipidomics Core Facility Main Medical Center Research Institute



The most consistent and flexible system available

Fully customizable infusion eliminates carry-over concerns while saving you time. Allows for sample infusion at flow rates between 5 and 100 µL/min - directly into your ESI source via the injection valve - precisely and without pulsation. Compatible with both shotgun and targeted lipid analysis.

The powerful Chronos software can execute important operational and integration functions such as reading and




Upgrade your sample throughput without worrying about carry-over

The PAL injector head is a syringe pump itself, making it the ideal sample delivery system for ESI flow rates, unlike conventional HPLC pumps.


Many lipids are notorious for their adherence to glass and other surfaces resulting in carryover between samples. The special LCMS syringe tool incorporates unique features to minimize this carry-over to below 1% for the abundant neutral and phospholipids and even lower for other species.



The illustration below shows how this rapid washing is achieved. After washing the inside, the system will wash the needle externally and be ready to aspirate the next sample.



Greatest experimental control with powerful, user-friendly software



An example Chronos infusion method

The infusion process is controlled using our Chronos software, allowing maximum compatibility with your preferred Chromatography Data System, LIMS and MS data systems.

Method parameters are available for flow rates, infusion times, stabilization delays, cleanup conditions, and MS acquisition method switching.

With an optional second valve, standard chromatography methods can be used in the same sequence with direct infusion samples.

High throughput, high sensitivity and low carry-over: the new infusion automation from Trajan LEAP Automation offers you everything you need to bring your lipidomics mass spec workflow to its peak performance.

Infusion Automation

Expert infusion automation for your high sensitivity lipidomics experiments - the specifications can be found in the table below.

Specifications	
Flow rate range	5 - 100 μL/min
Sample formats	Sealed and unsealed vials (1mL, 2mL, 10 mL); sealed and unsealed plates (96 and 384-well)
Optional features	Sample cooling drawers (4 - 40°C), inert flow path, 2nd valve for chromatography autosampling
Sample capacity	Expandable capacity, from three plate tray holder, to 6 or 12 plate drawer stack(s). Capacity depends on vial/ tray formats and drawer configuration, maximum 4 stacks.
Signal stability	Signal stability on a single sample has been demonstrated for 20 minutes. Long term stability confirmed in successful analysis of 400 samples over a 36-hour period (5 minute injection per sample).
Computer requirements (minimum)	Operating system: Windows 7 SP1 (32 or 64 bit), Windows 8.1 or Windows 10. Hardware: Intel Dual Core 2.0 GHz, 60 GB hard disk, 4 GB RAM, dedicated ethernet port or free location on local private network switch with 1 free USB port for security dongle. Monitor: 1366x768 (recommended 1440x900 or 1920x1080)

Visit us at <u>www.leaptec.com</u> or contact your regional Trajan LEAP Automation representative for assistance.

Trajan Scientific and Medical Science that benefits people. Trajan is actively engaged in developing and delivering solutions that have a positive impact on human wellbeing. Our vision revolves around collaborative partnerships that improve workflows, delivering better results.



Legal Statements

CTC Analytics AG reserves the right to make improvements and/or changes to the product(s) described in this document at any time without prior notice.

CTC Analytics AG makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.

Under no circumstances shall CTC Analytics AG be held liable for any coincidental damage or damages arising as a consequence of or from the use of this document.

© 2018 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

Imprint

Date of print: 12.2018

CTC Analytics AG Industriestrasse 20 CH-4222 Zwingen Switzerland T +41 61 765 81 00 Contact: info@ctc.ch