

NEW
for CTC PAL
R-Syringe

golden coloured plunger
end piece
indication for coated needle
protect the plunger before
damages

plunger guide
longer lifetime

Coloured
metal flange
Indication for
Needle I.D.


glassbody
DURAN®Glas
3.3 Borosilkat
polished

Plunger seal
PTFE black- seal
or PE-HD

no contact of the
liquids with the glue
because of the seal

needle
coated

Description	Seal	Part.No.
10µlR-SYR FN0,72(G22s)b51	PTFE – black-seal	2620440
25µlR-SYR FN0,72(G22s)b51	PTFE – black-seal	2620540
25µlR-SYR FN0,72(G22)b51	PTFE – black-seal	2620541
25µlR-SYR FN0,72(G22s)b51	PE-HD	2620542
25µlR- SYRFN0,72(G22)b51	PE-HD	2620543
50µlR-SYR FN0,72(G22s)b51	PTFE – black-seal	2620640
50µlR-SYR FN0,72(G22)b51	PTFE – black-seal	2620641
50µlR-SYR FN0,72(G22s)b51	PE-HD	2620642
50µlR- SYRFN0,72(G22)b51	PE-HD	2620643
100µlR-SYR FN0,72(G22s)b51	PTFE – black-seal	2620740
100µlR-SYR FN0,72(G22)b51	PTFE – black-seal	2620741
100µlR-SYR FN0,72(G22s)b51	PE-HD	2620742
100µlR- SYRFN0,72(G22)b51	PE-HD	2620743



G22
(blue)
I.D.
0,40mm

G22s
(red)
I.D.
0,20mm